



## HMX

---

Method number: PV2032

Target concentration: 0.2 mg/m<sup>3</sup> (arbitrary). There is no OSHA PEL or ACGIH TLV for HMX.

Procedure: Samples are collected by drawing air through glass fiber filters. Samples are extracted with acetone and analyzed by high performance liquid chromatography (HPLC) with a UV detector.

Recommended air volume and sampling rate: 500 minutes at 1 Lpm (500 L)

Detection limit of the analytical procedure: 10.5 ng/injection

Status of method: Partially Validated method. This method has been only partially evaluated and is presented for information and trial use.

July 1987 (final)

Yihlin Chan

Carcinogen and Pesticide Branch  
OSHA Salt Lake Technical Center  
Salt Lake City UT 84115-1802



2.1.2 Glass fiber filter, 37-mm diameter, Gelman Type A or equivalent.

2.1.3 Filter holder for 37-mm filters, Millipore M000037A0 or equivalent.

## 2.2 Sampling technique

Use standard air sampling procedure as specified in the OSHA Instruction CPL 2-2.20A, Chapter II: Standard Method for Sampling Air Contaminants.

## 2.3 Recommended air volume and sampling rate

2.3.1 The recommended air volume is 500 L.

2.3.2 The recommended sampling rate is 1 Lpm.

## 2.4 Extraction efficiency

Six glass fiber filters were each spiked with 23.44 ug of HMX. After overnight storage, two filters were extracted with 5 mL of acetone. The average recovery of HMX was 95.6%.

Table 2.4  
Extraction Efficiency

sample #	HMX recovered	% recovered
1	22.8	97.3
2	22.0	93.9

average = 95.6%

## 2.5 Retention efficiency

Humid air (90% RH, 500 L) was drawn through the remaining four filters in Section 2.4. Two filters were analyzed. The average recovery of HMX was 93.7%.

Table 2.5  
Retention Efficiency

sample #	HMX recovered	% recovered
3	21.4	96.0
4	20.5	91.3

average = 93.7%

## 2.6 Storage test (7 days)

Three glass fiber filters were each spiked with 23.44 ug of HMX. Humid air (RH 74%, 276 L at 1 Lpm) was pulled through the filters. The filters were stored at room temperature for 7 days, extracted, and analyzed. The average recovery of the HMX was 102.7%.

Table 2.6  
Storage Test (7-days)

sample #	HMX recovered	% recovered
1	24.21	103.3
2	24.15	103.0
3	23.86	101.8

average = 102.7%

## 2.7 Interferences (Sampling)

There are no known interferences to the sampling procedure.

## 3 Analytical procedure

### 3.1 Apparatus

3.1.1 High performance liquid chromatograph

3.1.2 Zorbax ODS HPLC column

3.1.3 UV detector. Waters 490 Programmable Multiwavelength Detector was used in this study.

3.1.4 Strip chart recorder

### 3.2 Reagents

3.2.1 HMX standard

3.2.2 Acetonitrile, HPLC grade

3.2.3 Water, HPLC grade

3.2.4 Acetone, reagent grade

### 3.3 Standard preparation

Weigh 3 to 5 mg of HMX in a 10-mL volumetric flask. Add acetone to the mark. Dilute standards to a working range of 0.7 to 20 ug/mL.

### 3.4 Sample preparation

Place glass fiber filter in a scintillation vial. Add 5 mL of acetone. Shake on a mechanical shaker for 30 minutes.

### 3.5 Analysis

#### 3.5.1 Instrument conditions

Column: Zorbax ODS column  
Eluent: 40% acetonitrile/60% water  
Flow rate: 1.0 mL/min  
Detector: 235 nm (primary), 225 nm  
Retention time: 11.1 min

3.5.2 Chromatogram (see Figure 1)

3.6 Interferences

3.6.1 Any compound that has the same retention time as the HMX and absorbs at 235 nm and 225 nm is interference.

3.6.2 Most interference can be circumvented by altering the chromatographic conditions.

3.6.3 Retention time alone is not proof of a chemical identity. Confirmation by other means should be sought when possible.

3.7 Calculations

3.7.1 A calibration curve is constructed by plotting concentration versus detector response.

3.7.2 The concentration of a sample is determined from the calibration curve.

3.7.3 The concentration of HMX is given by:

$$\text{mg} / \text{m}^3 = \frac{(\mu\text{g} / \text{mL}, \text{blank corrected})(5 \text{ mL})}{(\text{air volume}, \text{L})(\text{desorption efficiency}, \text{decimal})}$$

4 Recommendations for further study

The method should be fully validated.




Figure 1 HPLC Chromatogram of HMX




Figure 2 Calibration Curve of HMX

## 5 References

- 5.1 Registry of Toxic Effects of Chemical Substances 1983-84, DHHS (NIOSH) Publication No. 86-103.
- 5.2 Grayson, Martin, ed., Kirk-Othmer Encyclopedia of Chemical Technology, Volume 9, New York: John Wiley & Sons, 1980.